

Product datasheet

Anti-sRANKR antibody [R12-31] (Phycoerythrin) ab95654

1 图像

概述

产品名称	Anti-sRANKR抗体[R12-31] (Phycoerythrin)
描述	大鼠单克隆抗体[R12-31] to sRANKR (Phycoerythrin)
宿主	Rat
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Mouse sRANKR protein
阳性对照	Mouse sRANKR-transfected cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: Less than or equal to 0.09% Sodium Azide Constituents: PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
克隆编号	R12-31
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab95654** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
Flow Cyt		Use 0.5-1µg for 10 <sup>5-8</sup> cells. Staining the cell sample in a final volume of 100 µl is recommended.

## 靶标

**相关性** sRANKR is a receptor for TNFSF11/RANKL/TRANCE/OPGL; essential for RANKL-mediated osteoclastogenesis. It is involved in the regulation of interactions between T-cells and dendritic cells.

**细胞定位** Membrane; Single-pass type I membrane protein

## 图片



Flow cytometry staining of sRANKR-transfected cells with 0.5  $\mu$ g of Rat IgG2a ? Isotype Control PE (open histogram) or 0.5  $\mu$ g of ab95654 (filled histogram). Total viable cells were used for analysis.

Flow Cytometry - Anti-sRANKR antibody [R12-31]  
(Phycoerythrin) (ab95654)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors