

Product datasheet

Anti-ACSF3 antibody ab73352

1 图像

概述

产品名称	Anti-ACSF3抗体
描述	小鼠多克隆抗体to ACSF3
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human ACSF3 protein (AAH28399.1)
阳性对照	ACSF3 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Constituent: Whole serum
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab73352** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 64 kDa.

靶标

**相关性** Acyl-CoA (ACSF3) synthases catalyze the initial reaction in fatty acid metabolism, by forming a thioester with CoA and may have some preference toward very long chain substrates.

**细胞定位** Mitochondrial

## 图片



**All lanes :** Anti-ACSF3 antibody (ab73352) at 5 µg/ml

**Lane 1 :** ACSF3 transfected 293T cell lysate

**Lane 2 :** Non transfected cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** goat anti-mouse (H&L)-HRP conjugate at 1/2500 dilution

**Predicted band size:** 64 kDa

**Observed band size:** 64 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors