

Product datasheet

Anti-ALPK2 antibody ab111909

1 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-ALPK2抗体 |
| 描述 | 兔多克隆抗体to ALPK2 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Human |
| 免疫原 | A peptide from C terminal residues of Human ALPK2 protein (NM_052947.3). |
| 阳性对照 | Human fetal liver lysate. |

性能

| | |
|------|--|
| 形式 | Lyophilised:Reconstitute in 200ul Sterile Distilled Water. |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 1% BSA, 98.98% PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | ab111909 is purified by a peptide affinity column. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab111909** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|--|
| WB | | 1/500 - 1/1000. Detects a band of approximately 237 kDa (predicted molecular weight: 237 kDa). |
| IHC-P | | 1/100 - 1/500. |

靶标

| | |
|-------|---|
| 功能 | Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. |
| 序列相似性 | Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily. Contains 1 alpha-type protein kinase domain. Contains 2 Ig-like (immunoglobulin-like) domains. |

图片



Anti-ALPK2 antibody (ab111909) at 1/500 dilution + Human fetal liver lysate

Predicted band size: 237 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors