

## Product datasheet

# Anti-ATP10B antibody ab122216

### 1 图像

#### 概述

<b>产品名称</b>	Anti-ATP10B抗体
<b>描述</b>	兔多克隆抗体to ATP10B
<b>宿主</b>	Rabbit
<b>经测试应用</b>	<b>适用于:</b> IHC-P
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	antigen sequence, corresponding to amino acids 463-561 (RLETPKELDS DGEEWTQYQC LSF SARWAQD PATMRSQKGA QPLRRSQSAR VPIQGHYRQR SMGHRESSQP PVAFSSSIEK DVTPDKNLLT KVRDAALWL) of Human ATP10B. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	Human cerebral cortex tissue.

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab122216** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标

组织特异性	Found in brain and in low levels in testis.
序列相似性	Belongs to the cation transport ATPase (P-type) (TC 3.A.3) family. Type IV subfamily.
细胞定位	Membrane.

## 图片



ab122216 at 1/400 dilution staining ATP10B in Paraffin-embedded Human Cerebral cortex tissue by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors