

Product datasheet

Anti-AU1 tag antibody (Agarose) ab3846

概述

产品名称	Anti-AU1 tag抗体(Agarose)
描述	山羊多克隆抗体to AU1 tag (Agarose)
宿主	Goat
偶联物	Agarose
经测试应用	适用于: ELISA, WB, ICC
免疫原	Synthetic peptide: DTYRYI conjugated to KLH.
	Run BLAST with Run BLAST with
常规说明	Antibodies were coupled to agarose beads using a cyanogen bromide method.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 0.58% Sodium chloride
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab3846** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		
WB		

应用	Ab评论	说明
----	------	----

ICC

应用说明

ELISA: Use at a dilution of: Coating: 1/100 – 1/500
 Primary: 1:1000 – 1/30,000

ICC: 1/200 - 1/2000.

WB: Use at a dilution of: Chemiluminescent detection: 1/1,000 – 1/30,000

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors