

Product datasheet

Anti-Adracalin antibody ab96351

1 图像

概述

产品名称	Anti-Adracalin抗体
描述	兔多克隆抗体to Adracalin
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant protein fragment with a sequence corresponding to a region within amino acids 320 and 516 of Human Adracalin
阳性对照	HeLaS3 cell lines

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, PBS, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96351** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 60 kDa.

靶标

相关性 Adracalin (AAAS) is expressed in both neuroendocrine and cerebral structures and may function in the normal development of the peripheral and central nervous system. It localizes to nuclear pore complexes (NPCs), large multiprotein assemblies that are the sole sites of nucleocytoplasmic transport. Defects in AAAS are the cause of achalasia-addisonianism-alacrima syndrome (AAA syndrome); also known as triple-A syndrome or Allgrove syndrome.

细胞定位 nuclear pore

图片



Anti-Adracalin antibody (ab96351) at 1/1000 dilution + HeLa S3 whole cell lysate at 30 µg

Predicted band size: 60 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors