

Product datasheet

Anti-BCAS2 antibody [EPR4375] ab108330

重组 RabMAb

2 图像

概述

产品名称	Anti-BCAS2抗体[EPR4375]
描述	兔单克隆抗体[EPR4375] to BCAS2
宿主	Rabbit
经测试应用	适用于: WB, IP, ICC/IF 不适用于: Flow Cyt or IHC-P
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide within Human BCAS2 aa 1-100. The exact sequence is proprietary.
阳性对照	SH-SY5Y, MCF-7, and Neuro-2a cell lysates, HeLa cells
常规说明	

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#)

This product is a recombinant rabbit monoclonal antibody.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol, 0.05% BSA, 50% Tissue culture supernatant
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	EPR4375
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab108330** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000 - 1/10000. Predicted molecular weight: 26 kDa.
IP		1/10 - 1/100.
ICC/IF		1/250 - 1/500.
应用说明		Is unsuitable for Flow Cyt or IHC-P.

靶标


功能	Involved in mRNA splicing.
组织特异性	Ubiquitously expressed.
序列相似性	Belongs to the SPF27 family.
细胞定位	Nucleus > nucleolus.

图片



Immunofluorescent staining of HeLa cells using ab108330 at a dilution of 1/250.

Immunocytochemistry/ Immunofluorescence - Anti-BCAS2 antibody [EPR4375] (ab108330)



Western blot - Anti-BCAS2 antibody [EPR4375]

(ab108330)

All lanes : Anti-BCAS2 antibody [EPR4375]

(ab108330) at 1/1000 dilution

Lane 1 : SH-SY5Y cell lysate

Lane 2 : MCF7 cell lysate

Lane 3 : Neuro-2a cell lysate

Lysates/proteins at 10 µg per lane.

Predicted band size: 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors