

Product datasheet

Anti-BDH1 antibody ab103500

2 图像

概述

产品名称	Anti-BDH1抗体
描述	兔多克隆抗体to BDH1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-343 of Human BDH1 (NP_004042.1).
阳性对照	HeLa cell lysate and BDH1 transfected 293T cell line lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab103500** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 38 kDa.

靶标


序列相似性 Belongs to the short-chain dehydrogenases/reductases (SDR) family.

## 图片



Anti-BDH1 antibody (ab103500) at 1/500 dilution + HeLa cell lysate at 50  $\mu$ g

**Predicted band size:** 38 kDa



**All lanes :** Anti-BDH1 antibody (ab103500) at 1/500 dilution

**Lane 1 :** BDH1 transfected 293T cell line

**Lane 2 :** Non transfected 293T cell line

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 38 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors