

Product datasheet

Anti-BDP / ARID3B antibody - C-terminal ab76499

1 图像

概述

产品名称	Anti-BDP / ARID3B抗体- C-terminal
描述	兔多克隆抗体to BDP / ARID3B - C-terminal
宿主	Rabbit
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide selected from the C-terminal between 418-447 amino acid region of human BDP/ ARID3B conjugated to KLH
阳性对照	Human testis tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab76499 is purified by Saturated Ammonium Sulphate precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab76499** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/1000.

应用	Ab评论	说明
IHC-P		1/10 - 1/50.

## 靶标

功能	Transcription factor which may be involved in neuroblastoma growth and malignant transformation. Favors nuclear targeting of ARID3A.
组织特异性	Expressed in placenta, testis and leukocytes. Expressed in neuroblastoma. Present in K562 erythrocytic leukemia cells (at protein level).
序列相似性	Contains 1 ARID domain. Contains 1 REKLES domain.
细胞定位	Nucleus.

## 图片



Formalin-fixed and paraffin-embedded human testis tissue, stained with ab76499, at 1/10 dilution, using a peroxidase-conjugated secondary antibody and DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-BDP / ARID3B antibody - C-terminal (ab76499)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors