

Product datasheet

Anti-BRDG 1 antibody ab103646

3 图像

概述

产品名称	Anti-BRDG 1抗体
描述	兔多克隆抗体to BRDG 1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Full length Human BRDG 1 (NP_036240.1; amino acids 1-295).
阳性对照	Human kidney and Mouse intestine tissue lysates, BRDG1 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab103646** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 34 kDa.

靶标

功能 In BCR signaling, appears to function as a docking protein acting downstream of TEC and

participates in a positive feedback loop by increasing the activity of TEC.

**序列相似性**

Contains 1 PH domain.

Contains 1 SH2 domain.


**翻译后修饰**

Phosphorylated on tyrosine by TEC.

**细胞定位**


Cytoplasm.

**图片**




Anti-BRDG 1 antibody (ab103646) at 1  $\mu$ g/ml  
+ Human kidney tissue lysate at 50  $\mu$ g

**Predicted band size:** 34 kDa



Anti-BRDG 1 antibody (ab103646) at 1  $\mu$ g/ml  
+ Mouse intestine tissue lysate at 50  $\mu$ g

**Predicted band size:** 34 kDa



**All lanes :** Anti-BRDG 1 antibody (ab103646)  
at 1  $\mu$ g/ml

**Lane 1 :** BRDG 1 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors