

## Product datasheet

# Anti-Bovine Serum Albumin antibody [BSA-7G10] ab3781

★★★★☆ 3 Abreviews 3 References

### 概述

产品名称	Anti-Bovine Serum Albumin抗体[BSA-7G10]
描述	小鼠单克隆抗体[BSA-7G10] to Bovine Serum Albumin
宿主	Mouse
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Mouse, Rat, Cow 预测可用于: Rabbit, Human
免疫原	Full length native protein (purified) (Cow).

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 1.2114% Tris, 0.7507% Glycine, 2% Sucrose
纯度	Protein A purified
克隆	单克隆
克隆编号	BSA-7G10
骨髓瘤	unknown
同种型	IgG1
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab3781** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000 - 1/10000.

应用	Ab评论	说明
WB	★★★★☆	1/1200 - 1/2000.

## 靶标

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors