

Product datasheet

Anti-C6orf186 antibody ab122388

2 图像

概述

产品名称	Anti-C6orf186抗体
描述	兔多克隆抗体to C6orf186
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment : LEQIGQLIFE IHLHWPGFEV SGSDSSVVRF WYLLKELEQ KDFRLFHSYK DLSKPQLFLK KDIFNASSCY TLSWVNTRW, corresponding to C terminal amino acids 287-365 of Human C6orf186 (Q5JXM2) Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human duodenum tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab122388** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		1/1000 - 1/2500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

应用	Ab评论	说明
----	------	----

靶标


序列相似性	Belongs to the UPF0624 family.
细胞定位	Secreted.

图片



ab122388 at 1/1000 staining C6orf186 in formalin fixed paraffin embedded Human duodenum tissue using immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C6orf186 antibody (ab122388)



Immunohistochemical staining of Human duodenum shows strong cytoplasmic positivity with a granular pattern using ab122388

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C6orf186 antibody (ab122388)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors