

## Product datasheet

# Anti-CCNL2 antibody ab115547

### 1 图像

#### 概述

<b>产品名称</b>	Anti-CCNL2抗体
<b>描述</b>	兔多克隆抗体to CCNL2
<b>宿主</b>	Rabbit
<b>经测试应用</b>	<b>适用于:</b> WB, IP, ELISA, IHC-P, ICC/IF
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Rat, Chicken 
<b>免疫原</b>	Synthetic peptide: KAQARGLLPG GTQVLDGTSG FSPAPKLVES PKEGKGSKPS PLSVKNTKRR LEGAKKAKAD SPVNLPKGR ESRSR, corresponding to internal sequence amino acids 309-384 of Human CCNL2. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	Human small intestine tissue

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.27% Potassium chloride, 0.88% Sodium chloride
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用


Our [Abpromise guarantee](#) covers the use of **ab115547** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 58 kDa.
IP		Use at an assay dependent concentration.
ELISA		1/2000 - 1/10000.
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use at an assay dependent concentration.

## 靶标

## 图片



Immunohistochemical image of CCNL2 in formalin fixed, paraffin embedded Human small intestine tissue stained with ab115547 at 5 µg/ml, followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CCNL2 antibody (ab115547)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors