

Product datasheet

Anti-CD15 antibody [28] ab20137

3 References

概述

产品名称	Anti-CD15抗体[28]
描述	小鼠单克隆抗体[28] to CD15
宿主	Mouse
特异性	Specific to CD15. 28 is useful for leukaemia typing and purification of stem cells.
经测试应用	适用于: IP, Flow Cyt, IHC-Fr
种属反应性	与反应: Human
免疫原	Cell preparation: Monocytes, separated from other peripheral blood leucocytes (PBL) on fibronectin plates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	28
骨髓瘤	Sp2/0-Ag14
同种型	IgM
轻链类型	unknown

应用

Our [Abpromise guarantee](#) covers the use of **ab20137** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		

应用	Ab评论	说明
Flow Cyt		
IHC-Fr		
应用说明	<p>Flow Cyt: Use at an assay dependent dilution.</p> <p>IHC-Fr: Use at an assay dependent dilution.</p> <p>Note: samples should be acetone fixed.</p> <p>IP: Use at an assay dependent dilution.</p> <p>Not tested in other applications.</p> <p>Optimal dilutions/concentrations should be determined by the end user.</p>	
靶标		
相关性	<p>CD15 is a carbohydrate adhesion molecule (and not a protein) that mediates phagocytosis and chemotaxis. Synthesis is directed by FUT4 in lymphoid cells and mature granulocytes, and by FUT9 in promyelocytes and monocytes.</p>	
细胞定位	<p>Golgi Apparatus; Membrane-bound form in trans cisternae of Golgi.</p>	
Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"		

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors