

Product datasheet

Anti-CD2 antibody [TP1/31] (CF405M) ab117737

2 图像

概述

产品名称	Anti-CD2抗体[TP1/31] (CF405M)
描述	小鼠单克隆抗体[TP1/31] to CD2 (CF405M)
宿主	Mouse
偶联物	CF405M. Ex: 408nm, Em: 452nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	T cell blasts activated with PMA and ionomycin
阳性对照	Normal Human peripheral blood cells.
常规说明	CF405M (Abs/Em Max: 408/450nm). Direct replacement for Pacific Blue dye [®] , BD Horizon™ V450.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Note: Contains antibody stabilizer.
纯度	Immunogen affinity purified
克隆	单克隆
克隆编号	TP1/31
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab117737** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use 5µl for 10 ⁶ cells. ab126026 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标


功能	CD2 interacts with lymphocyte function-associated antigen (LFA-3) and CD48/BCM1 to mediate adhesion between T-cells and other cell types. CD2 is implicated in the triggering of T-cells, the cytoplasmic domain is implicated in the signaling function.
序列相似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

图片



Fluorescence profile of normal Human peripheral blood unstained (blue) or stained with ab117737 (yellow) at 5 µl/10⁶ cells.

Flow Cytometry - Anti-CD2 antibody [TP1/31]
(CF405M) (ab117737)



Fluorescence profile of normal Human peripheral blood unstained (blue histogram) or stained with ab117737 (yellow histogram) at 5 $\mu\text{l}/10^6$ cells.

Flow Cytometry - Anti-CD2 antibody [TP1/31]
(CF405M) (ab117737)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors