

Product datasheet

Anti-CD28 antibody [B-T3] (FITC) ab27310

1 References 1 图像

概述

产品名称	Anti-CD28抗体[B-T3] (FITC)
描述	小鼠单克隆抗体[B-T3] to CD28 (FITC)
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	Recognises the B7/BB1 receptor (B7/BB1 R)
经测试应用	<b>适用于:</b> Flow Cyt
种属反应性	<b>与反应:</b> Human
免疫原	U266 cell line
常规说明	Abcam is committed to meeting high standards of ethical manufacturing and has decided to discontinue this product by June 2019 as it has been generated by the ascites method. We are sorry for any inconvenience this may cause.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium azide Constituents: PBS, 5% BSA
纯度	Ion Exchange Chromatography
克隆	单克隆
克隆编号	B-T3
骨髓瘤	x63-Ag8.653
同种型	IgG2a


应用

Our [Abpromise guarantee](#) covers the use of **ab27310** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		
应用说明	Flow Cyt: Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood.	
	Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.	
靶标		
功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.	
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.	
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.	
细胞定位	Membrane.	

## 图片



A typical staining pattern with ab27310 of lymphocytes

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors