

Product datasheet

Anti-CD89 antibody [MIP7c] ab53357

1 References 1 图像

概述

产品名称	Anti-CD89抗体[MIP7c]
描述	小鼠单克隆抗体[MIP7c] to CD89
宿主	Mouse
特异性	ab53357 recognises the EC1 domain of the CD89 extracellular region.
经测试应用	适用于: Flow Cyt, Inhibition Assay
种属反应性	与反应: Human
免疫原	Soluble CD89

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	MIP7c
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab53357** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

Inhibition Assay

应用说明 Flow Cyt: Use 1µg/100µl as primary antibody.

Inhib: Use 4µg for 10⁹ cells.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	Binds to the Fc region of immunoglobulins alpha. Mediates several functions including cytokine production.
组织特异性	Isoform A.1, isoform A.2 and isoform A.3 are differentially expressed between blood and mucosal myeloid cells. Isoform A.1, isoform A.2 and isoform A.3 are expressed in monocytes. Isoform A.1 and isoform A.2 are expressed in alveolar macrophages; however only one isoform is expressed at alveolar macrophages surfaces.
序列相似性	Contains 2 Ig-like C2-type (immunoglobulin-like) domains.
细胞定位	Secreted and Cell membrane.

图片



Neutrophil phagocytosis of IgA immune complexes in presence or absence (control) of ab53357.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors