

Product datasheet

Anti-CNOT3 antibody [2192C2a] ab51335

概述

产品名称	Anti-CNOT3抗体[2192C2a]
描述	小鼠单克隆抗体[2192C2a] to CNOT3
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Human
免疫原	Recombinant human CNOT3 fragment

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, 8mM Sodium phosphate, 3mM Potassium chloride, 140mM Sodium chloride, 1.5mM Potassium phosphate, pH 7.4
纯度	Protein G purified
纯化说明	ab51335 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.
克隆	单克隆
克隆编号	2192C2a
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab51335** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 82 kDa.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution.

靶标

功能	The CCR4-NOT complex functions as general transcription regulation complex.
组织特异性	Ubiquitous. Highly expressed in brain, heart, thymus, spleen, kidney, liver, small intestine, lung and peripheral blood leukocytes.
序列相似性	Belongs to the CNOT2/3/5 family.
细胞定位	Cytoplasm. Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors