

Product datasheet

CNPase overexpression 293T lysate (whole cell) ab94074

2 图像

概述

<b>产品名称</b>	CNPase overexpression 293T裂解物(whole cell)
<b>常规说明</b>	ab94074 is a 293T cell transfected lysate in which Human CNPase has been transiently over-expressed using a pCMV-CNPase plasmid. The lysate is provided in 1X Sample Buffer. <b>Note:</b> For more details about how the transfected lysate was prepared <a href="#">view preparation notes</a>
<b>经测试应用</b>	<b>适用于:</b> WB

性能


<b>形式</b>	Liquid
<b>存放说明</b>	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	pH: 6.80 Constituent: 100% 1x Sample Buffer
<b>纯度</b>	Whole Cell Lysate
<b>背景</b>	Similarity: Belongs to the cyclic nucleotide phosphodiesterase family.

应用


Our [Abpromise guarantee](#) covers the use of **ab94074** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



ab94074 at 15µg/lane on an SDS-PAGE gel.



**All lanes :** Anti-CNPase antibody (ab92914)  
at 1/500 dilution

**Lane 1 :** CNPase overexpression 293T lysate  
(whole cell) (ab94074)

**Lane 2 :** 293T non transfected lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L)  
HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors