

Product datasheet

Anti-COX1 / Cyclooxygenase 1 antibody ab52037

1 图像

概述

产品名称	Anti-COX1 / Cyclooxygenase 1抗体
描述	兔多克隆抗体to COX1 / Cyclooxygenase 1
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide within Human COX1/ Cyclooxygenase 1 aa 550 to the C-terminus (C terminal). The exact sequence is proprietary. Database link: P23219
阳性对照	Esophagus carcinoma

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 1% BSA, 0.0268% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab52037** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab 评论	说明

应用	Ab 评论	说明
IHC-P		1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. 1/100 for 30 min at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min.

靶标

功能	May play an important role in regulating or promoting cell proliferation in some normal and neoplastically transformed cells.
通路	Lipid metabolism; prostaglandin biosynthesis.
序列相似性	Belongs to the prostaglandin G/H synthase family. Contains 1 EGF-like domain.
细胞定位	Microsome membrane. Endoplasmic reticulum membrane.

图片



Human esophagus carcinoma stained with Anti-COX1, ab52037 at 1/100 for 30 min at RT

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-COX1 / Cyclooxygenase 1 antibody (ab52037)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors