

Product datasheet

Anti-Carboxypeptidase A antibody ab63806

1 图像

概述

产品名称	Anti-Carboxypeptidase A抗体
描述	兔多克隆抗体to Carboxypeptidase A
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat, Human, Pig
免疫原	Synthetic peptide surrounding amino acid 30 of human Carboxypeptidase A
阳性对照	Jurkat cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 30% Glycerol, 0.5% BSA, PBS, pH 7.2
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab63806** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 47 kDa (predicted molecular weight: 47 kDa).

靶标

功能	Carboxypeptidase that catalyzes the release of a C-terminal amino acid, but has little or no action with -Asp, -Glu, -Arg, -Lys or -Pro.
序列相似性	Belongs to the peptidase M14 family.
细胞定位	Secreted.

图片



Western blot - Anti-Carboxypeptidase A antibody (ab63806)

All lanes : Anti-Carboxypeptidase A antibody (ab63806) at 4 µg/ml

All lanes : Jurkat cell lysate

Predicted band size: 47 kDa

Observed band size: 47 kDa

Additional bands at: 30 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors