

Product datasheet

Cytochrome P450 2E1 overexpression 293T lysate (whole cell)
ab94254

2 图像

概述

产品名称	Cytochrome P450 2E1 overexpression 293T裂解物(whole cell)
常规说明	ab94254 is a 293T cell transfected lysate in which Human Cytochrome P450 2E1 has been transiently over-expressed using a pCMV-Cytochrome P450 2E1 plasmid. The lysate is provided in 1X Sample Buffer. Note: For more detailed how the transfected lysate was prepared view preparation notes
经测试应用	适用于: WB

性能

形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 6.80 Constituent: 100% 1x Sample Buffer
纯度	Whole Cell Lysate
背景	Function: Metabolizes several precarcinogens, drugs, and solvents to reactive metabolites. Inactivates a number of drugs and xenobiotics and also bioactivates many xenobiotic substrates to their hepatotoxic or carcinogenic forms. Similarity: Belongs to the cytochrome P450 family.


应用

Our [Abpromise guarantee](#) covers the use of **ab94254** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



ab94254 at 15µg/lane on an SDS-PAGE gel



All lanes : Anti-Cytochrome P450 2E1 antibody (ab90564) at 1/500 dilution

Lane 1 : Cytochrome P450 2E1 overexpression 293T lysate (whole cell) (ab94254)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors