

Product datasheet

Anti-DDDDK tag antibody [M2] - (Equivalent to FLAG antibodies from Sigma) (Alkaline Phosphatase) ab81655

概述

产品名称	Anti-DDDDK tag抗体[M2] - (Equivalent to FLAG抗体from Sigma) (Alkaline Phosphatase)
描述	小鼠单克隆抗体[M2] to DDDDK tag - (Equivalent to FLAG抗体from Sigma) (Alkaline Phosphatase)
宿主	Mouse
偶联物	Alkaline Phosphatase
特异性	ab81655 recognizes the DYKDDDDK tag sequence at an N-terminal, Met-N-terminal, C-terminal, or internal site.
经测试应用	适用于: WB
种属反应性	N/A
免疫原	Synthetic peptide corresponding to DDDDK tag . Sequence: DYKDDDDK

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.013% Sodium azide Constituents: 0.3025% Tris, 0.87% Sodium chloride
克隆	单克隆
克隆编号	M2
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab81655** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/5000.

靶标

相关性 This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors