

Product datasheet

Anti-DGKK antibody ab111042

2 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-DGKK抗体  |
| 描述    | 兔多克隆抗体to DGKK  |
| 宿主    | Rabbit   |
| 经测试应用 | 适用于: IHC-P, ICC/IF   |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Synthetic peptide derived from the C-terminal of Human DGKK. |
| 阳性对照  | COS7 cells; Human brain tissue.                              |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.  |
| 存储溶液 | pH: 7.40<br>Preservative: 0.02% Sodium azide<br>Constituents: 49.1% PBS, 0.88% Sodium chloride, 50% Glycerol<br>Note: without Mg <sup>2+</sup> and Ca <sup>2+</sup> |
| 纯度   | Immunogen affinity purified   |
| 纯化说明 | ab111042 was affinity-purified from Rabbit antiserum by affinity-chromatography using epitope-specific immunogen.   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab111042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|----|
|    |      |    |


| 应用     | Ab评论 | 说明   |
|--------|------|--|
| IHC-P  |      | 1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |
| ICC/IF |      | 1/100 - 1/500.   |

## 靶标


**相关性** DGKK (Diacylglycerol kinase kappa) phosphorylates diacylglycerol (DAG) to generate phosphatidic acid (PA).

**细胞定位** Cytoplasm. Membrane; Peripheral membrane protein

## 图片



ab111042 at 1/100 dilution staining DGKK in COS7 cells by Immunofluorescence. The image on the right is treated with the synthesized peptide.



ab111042 at 1/50 dilution staining DGKK in paraffin-embedded Human brain tissue by Immunohistochemistry. The image on the right is treated with the synthesized peptide.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors