

Product datasheet

Anti-Hsp40 antibody ab80160

1 References

概述

产品名称	Anti-Hsp40抗体
描述	兔多克隆抗体to Hsp40
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Chlamydomonas reinhardtii
免疫原	Recombinant full length protein corresponding to Chlamydomonas reinhardtii Hsp40.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab80160** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/5000. Predicted molecular weight: 41 kDa.

靶标

功能	Interacts with HSP70 and can stimulate its ATPase activity. Stimulates the association between HSC70 and HIP.
序列相似性	Contains 1 J domain.

细胞定位

Cytoplasm. Nucleus. Nucleus > nucleolus. Translocates rapidly from the cytoplasm to the nucleus, and especially to the nucleoli, upon heat shock.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors