

## Product datasheet

# Donkey F(ab')<sub>2</sub> Anti-Rabbit IgG - H&L (DyLight® 594), pre-adsorbed ab98490

1 References 1 图像

### 概述

产品名称	驴F(ab') <sub>2</sub> 抗兔IgG - H&L (DyLight® 594), 预吸附二抗	
宿主	Donkey	
靶标种属	Rabbit	
经测试应用	适用于: IHC-P, ICC/IF, Flow Cyt	
最小交叉反应	Chicken, Cow, Goat, Horse, Human, Mouse, Rat, Sheep	<a href="#">more details</a>
偶联物	DyLight® 594. Ex: 593nm, Em: 618nm	

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 594.
克隆	多克隆
同种型	IgG

### 应用


Our [Abpromise guarantee](#) covers the use of **ab98490** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/50 - 1/500.

应用	Ab评论	说明
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

## 图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

DyLight®-Donkey F(ab')<sub>2</sub> polyclonal Secondary Antibody to Rabbit IgG - H&L (DyLight® 594), pre-adsorbed(ab98490)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors