




Product datasheet

Anti-UBE2O antibody ab118832

1 图像

概述

<b>产品名称</b>	Anti-UBE2O抗体
<b>描述</b>	山羊多克隆抗体to UBE2O
<b>宿主</b>	Goat
<b>经测试应用</b>	<b>适用于:</b> ELISA, IHC-P
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Cow, Dog, Pig 
<b>免疫原</b>	Synthetic peptide corresponding to Human UBE2O aa 1279-1292 (C terminal). Sequence: ALLEAGMPECTEDK  Database link: <a href="#">NP_071349.2</a>   <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
<b>阳性对照</b>	Human skeletal muscle tissue.
<b>常规说明</b>	Previously labelled as E2 230K.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 99% Tris buffered saline, 0.5% BSA
<b>纯度</b>	Protein G purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab118832** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/32000.
IHC-P		Use a concentration of 3.75 µg/ml.

## 靶标

功能	Catalyzes the covalent attachment of ubiquitin to other proteins.
组织特异性	Predominantly expressed in skeletal muscle and heart.
通路	Protein modification; protein ubiquitination.
序列相似性	Belongs to the ubiquitin-conjugating enzyme family.

## 图片



ab118832, at 3.75 µg/ml, staining UBE2O in formalin fixed, paraffin embedded Human skeletal muscle tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-UBE2O antibody (ab118832)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors