

Product datasheet

Anti-EN1 / Engrailed 1 antibody [EPR3282] ab108598

重组 RabMAb

1 图像

概述

产品名称	Anti-EN1 / Engrailed 1 抗体[EPR3282]
描述	兔单克隆抗体[EPR3282] to EN1 / Engrailed 1
宿主	Rabbit
经测试应用	适用于: WB, IP, Flow Cyt 不适用于: ICC or IHC-P
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide within Human EN1/ Engrailed 1. The exact sequence is proprietary.
阳性对照	U87-MG, SH-SY5Y, 293T and Human fetal brain lysates
常规说明	

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#)

This product is a recombinant rabbit monoclonal antibody.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol, 9.85% Tris glycine, 50% Tissue culture supernatant
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	EPR3282
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab108598** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000 - 1/10000. Detects a band of approximately 55 kDa (predicted molecular weight: 40 kDa).
IP		1/10 - 1/100.
Flow Cyt		1/100 - 1/500. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.


应用说明 Is unsuitable for ICC or IHC-P.

靶标

序列相似性 Belongs to the engrailed homeobox family.
Contains 1 homeobox DNA-binding domain.

细胞定位 Nucleus.

图片



All lanes : Anti-EN1 / Engrailed 1 antibody
[EPR3282] (ab108598) at 1/1000 dilution

Lane 1 : U87-MG lysate

Lane 2 : SH-SY5Y lysate

Lane 3 : 293T lysate

Lane 4 : Human fetal brain lysate

Lysates/proteins at 10 µg per lane.

Predicted band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors