

Product datasheet

Anti-ETL antibody ab12989

1 图像

概述

产品名称	Anti-ETL抗体
描述	兔多克隆抗体to ETL
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: PBS, pH 7.7
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab12989** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.

靶标

功能	Could be involved in cardiac development.
----	---

<b>组织特异性</b>	Mainly expressed in smooth muscle.
<b>序列相似性</b>	Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily. Contains 2 EGF-like domains. Contains 1 GPS domain.
<b>发展阶段</b>	Up-regulated in the adult heart.
<b>结构域</b>	The transmembrane domain is not required for cleavage, but it is required for dimer formation.
<b>翻译后修饰</b>	Proteolytically cleaved into 2 subunits, an extracellular alpha subunit and a seven-transmembrane subunit.
<b>细胞定位</b>	Cell membrane.

## 图片



Ab12989 positively staining ETL in formalin fixed, paraffin embedded heart tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ETL antibody (ab12989)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors