

Product datasheet

# Anti-Ethinylestradiol antibody ab59670

概述

|       |   |
|-------|---|
| 产品名称  | Anti-Ethinylestradiol抗体                             |
| 描述    | 羊多克隆抗体to Ethinylestradiol                           |
| 宿主    | Sheep   |
| 经测试应用 | 适用于: ELISA  |
| 种属反应性 | Not applicable                                      |
| 免疫原   | Ethinylestradiol conjugated to Bovine thyroglobulin |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.     |
| 存储溶液 | Preservative: 0.09% Sodium Azide<br>Constituents: 20mM PBS, 150mM Sodium chloride, pH 7.2 |
| 纯度   | IgG fraction  |
| 纯化说明 | Ig fraction prepared by caprylic acid and ammonium sulphate precipitation                 |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab59670** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用    | Ab评论 | 说明                                 |
|-------|------|------------------------------------|
| ELISA |      | Use a concentration of 1.25 µg/ml. |

靶标

相关性 Ethinylestradiol, (8R,9S,13S,14S,17R)-17-ethynyl-13-methyl-7,8,9,11,12,14,15,16-octahydro-6H-cyclopenta[a]phenanthrene-3,17-diol, is a semisynthetic alkylated Estradiol with a 17-alpha-

ethinyl substitution. It has high estrogenic potency when administered orally, and is often used as the estrogenic component in oral contraceptives. It has a molecular weight 296.40336 g/mol and molecular formula  $C_{20}H_{24}O_2$ .

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors