

Product datasheet

Anti-FPGT antibody ab97914

2 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-FPGT抗体 |
| 描述 | 兔多克隆抗体to FPGT |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant protein fragment, corresponding to a region within amino acids 146-423 of Human FPGT (NP_003829). |
| 阳性对照 | 293T, A431, H1299, HeLaS3 or Raji cell lines. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 1.21% Tris, 0.75% Glycine |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab97914** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| WB | | 1/500 - 1/3000. Predicted molecular weight: 67 kDa. |
| IHC-P | | 1/100. |

靶标

| | |
|-------|---|
| 功能 | Catalyzes the formation of GDP-L-fucose from GTP and L-fucose-1-phosphate. Functions as a salvage pathway to reutilize L-fucose arising from the turnover of glycoproteins and glycolipids. |
| 组织特异性 | Expressed in many tissues. |
| 细胞定位 | Cytoplasm. |


图片



Anti-FPGT antibody (ab97914) at 1/500 dilution + Raji whole cell lysate at 30 μ g

Predicted band size: 67 kDa

7.5% SDS PAGE



ab97914, at 1/100 dilution, staining FPGT in MDA-MB-468 xenograft tissue by Immunohistochemistry. Formalin-fixed, paraffin-embedded section.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors