

Product datasheet

Anti-GIP antibody ab22624

★★★★☆ 1 Abreviews

概述

产品名称	Anti-GIP抗体
描述	兔多克隆抗体to GIP
宿主	Rabbit
特异性	ab22624 reacts with Glucose-dependent Insulinotropic Peptide (GIP).
经测试应用	适用于: IHC-P, IHC-R, Flow Cyt
种属反应性	与反应: Mouse, Pig
免疫原	Full length native protein (purified) (Pig) conjugated to KLH by a CDI linker.
阳性对照	Mammalian duodenal (K) endocrine cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab22624** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.

应用	Ab评论	说明
IHC-R		1/1000 - 1/8000. 1/1000 - 1/8000 (using method A) and 1/200-1/1200 (using method B) Dilute to working strength with phosphatebuffered saline (PBS) containing 0.1% w/v sodium azide. Method A: using overnight incubation and indirect immunoperoxidase procedures, such as PAP and ABCperoxidase. Method B: using overnight incubation and indirect immunofluorescence procedures.
Flow Cyt	★★★★★	1/150. ab171870 - Rabbit polyclonal IgG, is suitable for use as an isotype control with this antibody.

靶标

功能	Potent stimulator of insulin secretion and relatively poor inhibitor of gastric acid secretion.
序列相似性	Belongs to the glucagon family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors