

Product datasheet

Anti-GLIS1 antibody ab96838

1 图像

概述

产品名称	Anti-GLIS1 抗体
描述	兔多克隆抗体 to GLIS1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide containing a sequence corresponding to a region within amino acids 559-620 (HSLSTPLPAT GYEALAEASC PTALPQQPSE DVVSSGPEDC GFFPNGAFDH CLGHIPSIYT DT) of Human GLIS1 (NP_671726). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Jurkat whole cell lysate; Hela cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 2.68% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96838** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 66 kDa.

## 靶标

功能	Acts as both a repressor and activator of transcription. Binds to the consensus sequence 5'-GACCACCCAC-3'.
序列相似性	Belongs to the GLI C2H2-type zinc-finger protein family. Contains 5 C2H2-type zinc fingers.
细胞定位	Nucleus.
形式	There are 2 isoforms produced by alternative splicing.

## 图片



Anti-GLIS1 antibody (ab96838) at 1/1000 dilution + Jurkat whole cell lysate at 30 µg

**Predicted band size:** 66 kDa

7.5% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors