

## Product datasheet

# GOLGA6 peptide (536-547) ab45812

### 概述

**产品名称** GOLGA6多肽(536-547)

### 描述

**性质** Synthetic

### 氨基酸序列

**种属** Human

**序列** C-DLPKEKADGTEQ

**氨基酸** 536 to 547

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab45812** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**应用** Blocking

**形式** Liquid

### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### 常规信息

**组织特异性** Highly expressed in seminiferous tubes in testis. Highly expressed in spermatids, barely detectable in late pachytene spermatocytes, and not detectable in spermatogonia. Detected at intermediate levels in pancreas and lymph nodes, and at much lower levels in spleen, peripheral blood leukocytes, skeletal muscle, liver, lung, placenta, brain and heart.

**序列相似性** Belongs to the GOLGA6 family.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors