

## Product datasheet

# Goat Anti-Armenian hamster IgG H&L ab5738

### 概述

产品名称	山羊抗亚美尼亚仓鼠IgG H&L
宿主	Goat
靶标种属	Armenian hamster
经测试应用	适用于: ELISA, WB, ICC/IF, IP, Immunodiffusion, Conjugation
免疫原	Armenian Hamster IgG whole molecule.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s). Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Armenian Hamster IgG and Armenian Hamster Serum.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab5738** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/20000 - 1/100000.
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.

应用	Ab评论	说明
IP		Use at an assay dependent dilution.
Immunodiffusion		Use at an assay dependent dilution.
Conjugation		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors