

Product datasheet

Goat Anti-Chicken IgY H&L (10nm Gold) ab41511

1 References

概述

| | |
|-------|---------------------------------|
| 产品名称 | 山羊抗鸡IgY H&L (10nm Gold) |
| 宿主 | Goat |
| 靶标种属 | Chicken |
| 经测试应用 | 适用于: Electron Microscopy |
| 免疫原 | Chicken IgG whole molecule. |
| 偶联物 | Gold 10nm |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | Preservative: 15mM Sodium Azide Constituents: 1% BSA, PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system [Moeremans et al., J. Immunol., Methods 74, (1984), 353.] |
| 克隆 | 多克隆 |
| 同种型 | IgG |
| 常规说明 | This product is an immunogold reagent built around colloidal gold particles of 10nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 10nm particle size conjugate shows less than 10% variation. |

应用

Our [Abpromise guarantee](#) covers the use of **ab41511** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|---------------------|------|----|
| Electron Microscopy | | |

应用说明

Electron Microscopy: Use at an assay dependent dilution. The method of production of this product allows the small particle size of 10nm to be used instead of 5nm with the same sensitivity but increased visibility in the standard transmission electron microscope (See also ven Bergen en Henegouwen and Leunissen, Histochemistry, 85,(1986), 81.)

Light Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors