

Product datasheet

Goat F(ab')₂ Anti-Rat IgG - (Fab')₂ (Biotin), pre-adsorbed ab6275

1 References

概述

产品名称	山羊F(ab') ₂ 抗大鼠IgG - (Fab') ₂ (Biotin),预吸附二抗
宿主	Goat
靶标种属	Rat
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Biotin, anti-Goat Serum, Rat IgG, Rat IgG F(ab') ₂ and Rat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Rat IgG F(c) or Bovine, Horse and Human Serum Proteins.
经测试应用	适用于: ICC/IF, Dot blot, WB, ELISA, IHC, Immunomicroscopy
最小交叉反应	Cow, Horse, Human more details
免疫原	Rat IgG F(ab') ₂ fragment
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 1.46% Sodium chloride, 1.5% BSA, 0.164% Sodium phosphate BSA is IgG and protease free
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Goat IgG molecule
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab6275** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/500 - 1/5000.
Dot blot		Use at an assay dependent dilution.
WB		1/2000 - 1/10000.
ELISA		1/20000 - 1/400000.
IHC		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors