

Product datasheet

Goat Anti-Human IgG H&L (HRP) ab6858

9 References 1 图像

概述

产品名称	山羊抗人IgG H&L (HRP)
描述	山羊多克隆抗体二抗到人IgG - H&L (HRP)
宿主	Goat
靶标种属	Human
经测试应用	适用于: ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB
免疫原	human IgG whole molecule
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Horseradish Peroxidase (HRP)
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab6858** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.

应用	Ab评论	说明
Dot blot		Use at an assay dependent concentration.
ELISA		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.
WB		1/5000 - 1/25000.

## 图片



Dot blot analysis of human IgG using ab6858 at 1/1000 dilution.

The indicated quantities of human IgG were spotted. Blocking was performed for 1 hour at room temperature using a proprietary blocking buffer prior to addition of ab6858.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors