

Product datasheet

Goat Anti-Human IgG H&L (TRITC) ab6855

概述

| | |
|-------|--|
| 产品名称 | 山羊抗人IgG H&L (TRITC) |
| 宿主 | Goat |
| 靶标种属 | Human |
| 经测试应用 | 适用于: Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA |
| 免疫原 | human IgG whole molecule |
| 偶联物 | TRITC. Ex: 547nm, Em: 572nm |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 7.60 Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Potassium phosphate |
| 纯度 | Affinity purified |
| 纯化说明 | This product was prepared from monospecific antiserum by immunoaffinity chromatography using human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. |
| 共轭说明 | Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.8 moles TRITC per mole of Rabbit IgG |
| 克隆 | 多克隆 |

应用

Our [Abpromise guarantee](#) covers the use of **ab6855** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|------------------|------|----|
| Immunomicroscopy | | |

| 应用 | Ab评论 | 说明 |
|-------------|---|----|
| Flow Cyt | | |
| ICC/IF | | |
| IHC-P | | |
| IHC-Fr | | |
| ELISA | | |
| 应用说明 | Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency | |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors