

## Product datasheet

# Goat Anti-Mouse IgG2a (SureLight® Allophycocyanin) ab72558

### 概述

产品名称	山羊抗小鼠IgG2a (SureLight® Allophycocyanin荧光)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
免疫原	Mouse IgG2a
偶联物	SureLight® Allophycocyanin. Ex: 652nm, Em: 657nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
纯度	Size exclusion
纯化说明	Purified using size exclusion chromatography as well as some proprietary methods for purifying the final materials.
共轭说明	Excitation max. $\lambda$ 650 nm. Emission max. $\lambda$ 660 nm.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab72558** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors