

Product datasheet

Goat Anti-Mouse IgM mu chain (Alkaline Phosphatase) preadsorbed ab98672

概述

产品名称	山羊抗小鼠IgM mu chain (Alkaline Phosphatase)预吸附二抗	
宿主	Goat	
靶标种属	Mouse	
经测试应用	适用于: IHC-P, ELISA, WB, ICC	
最小交叉反应	Human, Rat	more details
偶联物	Alkaline Phosphatase	

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, 0.1M Sodium chloride, 50mM HEPES, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 7.1
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab98672** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/100 - 1/1000.
ELISA		1/1000 - 1/10000. as primary

应用	Ab评论	说明
WB		1/1000 - 1/10000. Chemiluminescent detection.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors