

Product datasheet

Anti-HODE antibody ab28204

概述

产品名称	Anti-HODE抗体
描述	山羊多克隆抗体to HODE
宿主	Goat
特异性	REACTIVITY: HAPTEN REACTIVITY 9(±)-HODE 100.0% 9(S)-HODE 100.0% 9(R)-HODE 100.0% 13(S)-HODE 1.2% 13(R)-HODE 1.2% 9-oxo-ODE 1.2% 13-oxo-ODE 2.4% 11(S)-HETE 0.0% 15(S)-HETE 0.0% Linoleic Acid 0.0%
经测试应用	适用于: RIA, ELISA
免疫原	Chemical/ Small Molecule: 9(S)-HODE; 9(R)-HODE.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: PBS
纯化说明	Ammonium sulphate fractionation.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab28204** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
RIA		
ELISA		

应用说明	ELISA: Use at a concentration of 2 µg/well. RIA: Use at a concentration of 2 µg/well.
------	------------------------------------------------------------------------------------------

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

The hydroxyoctadecadienoic acid (HODE) derivatives of linoleic acid, 9(R)-HODE, 9(S)-HODE and 13(S)-HODE, are the most widely distributed of the known linoleic acid metabolites. These compounds exhibit interesting biologic activities, including the regulation of platelet function, maintenance of vascular thromboresistance and transduction of the cellular responses to certain growth factors. HODE derivatives may also influence certain pathological states including, psoriasis, the development of atherosclerosis and the development of cancer. 13-S-Hydroxyoctadecadienoic acid (13-S-HODE) is the product of 15 lipoxygenase (15 LOX) metabolism of linoleic acid. It is thought to enhance cellular mitogenic response to some growth factors.

细胞定位

Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors