

Product datasheet

Anti-Hemopexin antibody (Biotin) ab51515

概述

产品名称	Anti-Hemopexin抗体(Biotin)
描述	兔多克隆抗体to Hemopexin (Biotin)
宿主	Rabbit
偶联物	Biotin
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Human Hemopexin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 50% Glycerol, PBS, pH 7.5
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab51515** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use a concentration of 4 µg/ml.

靶标

功能 Binds heme and transports it to the liver for breakdown and iron recovery, after which the free

	hemopexin returns to the circulation.
组织特异性	Expressed by the liver and secreted in plasma.
序列相似性	Belongs to the hemopexin family. Contains 5 hemopexin-like domains.
翻译后修饰	N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors