


Product datasheet

Anti-Hepsin antibody ab27765

概述

产品名称	Anti-Hepsin抗体
描述	兔多克隆抗体to Hepsin
宿主	Rabbit
特异性	ab27765 recognises hepsin.
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Synthetic peptide: EWIFQAIKTHSEASGMVTQL , corresponding to C terminal amino acids 398-417 of Human Hepsin Run BLAST with Run BLAST with
阳性对照	Recombinant Human Hepsin protein (ab114907) can be used as a positive control in WB.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 1% BSA
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab27765** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/40.

靶标

功能	Plays an essential role in cell growth and maintenance of cell morphology.
组织特异性	Present in most tissues, with the highest level in liver.
序列相似性	Belongs to the peptidase S1 family. Contains 1 peptidase S1 domain. Contains 1 SRCR domain.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors