

Product datasheet

Human DUSP8 peptide ab22818

概述

产品名称 人DUSP8多肽

描述

性质 Synthetic

氨基酸序列

种属 Human

序列 AGDRLPRKVMDAK-C

氨基酸 2 to 14

技术指标

Our [Abpromise guarantee](#) covers the use of **ab22818** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking - Blocking peptide for Anti-DUSP8 antibody ([ab1363](#))

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

常规信息

功能 This protein shows both activity toward tyrosine-protein phosphate as well as with serine/threonine-protein phosphate.

组织特异性 Abundant in brain, heart and skeletal muscle.

序列相似性 Belongs to the protein-tyrosine phosphatase family. Non-receptor class dual specificity subfamily.

Contains 1 rhodanese domain.

Contains 1 tyrosine-protein phosphatase domain.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors