

## Product datasheet

# Human FSTL1 peptide ab23213

### 1 References

#### 概述

**产品名称** 人FSTL1多肽

#### 描述

**性质** Synthetic

**氨基酸序列**

**种属** Human

**序列** C-TAEKTKRVSTKEI

**氨基酸** 296 to 308

#### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab23213** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**应用** Blocking - Blocking peptide for Anti-FSTL1 antibody ([ab11805](#))

**形式** Liquid

#### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

#### 常规信息

**功能** May modulate the action of some growth factors on cell proliferation and differentiation. Binds heparin.

**组织特异性** Overexpressed in synovial tissues from rheumatoid arthritis.

**序列相似性** Contains 2 EF-hand domains.  
Contains 1 follistatin-like domain.  
Contains 1 Kazal-like domain.  
Contains 1 VWFC domain.

**细胞定位** Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors