

## Product datasheet

# Human PKR peptide ab51484

### 概述

产品名称 人PKR多肽

### 描述

性质 Synthetic

### 氨基酸序列

种属 Human

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab51484** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

### 制备和贮存

稳定性和存储 Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

pH: 7.20

Preservative: 0.02% Sodium azide

Constituents: PBS, 0.754% Sodium chloride, 0.1% BSA

### 常规信息

功能 Following activation by double-stranded RNA in the presence of ATP, the kinase becomes autophosphorylated and can catalyze the phosphorylation of the translation initiation factor EIF2S1, which leads to an inhibition of the initiation of protein synthesis. Double-stranded RNA is generated during the course of a viral infection.

序列相似性 Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. GCN2 subfamily. Contains 2 DRBM (double-stranded RNA-binding) domains. Contains 1 protein kinase domain.

翻译后修饰 Autophosphorylated on several Ser and Thr residues. Autophosphorylation of Thr-451 is

dependent on Thr-446 and is stimulated by dsRNA binding and dimerization.  
Autophosphorylation apparently leads to the activation of the kinase.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors