


Product datasheet

Anti-IL12 antibody ab124635

1 References 1 图像

概述

产品名称	Anti-IL12抗体
描述	山羊多克隆抗体to IL12
宿主	Goat
特异性	ab124635 was raised against a synthetic peptide from human IL12A / p35
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human 预测可用于: Dog, Pig 
免疫原	Synthetic peptide corresponding to residues in Human IL12.
阳性对照	Human intestine (Paneth cells) tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 99% Tris buffered saline, 0.5% BSA
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab124635** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/16000. Optimal dilution to be determined by the user.

应用	Ab评论	说明
IHC-P		Use a concentration of 3.75 µg/ml.

靶标

功能	Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.
序列相似性	Belongs to the IL-6 superfamily.
细胞定位	Secreted.

图片



ab124635, at 3.75µg/ml, staining IL12 in formalin-fixed, paraffin-embedded Human intestine (paneth cells) tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-IL12 antibody (ab124635)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors