

Product datasheet

Anti-ILDR2 antibody ab121369

3 图像

概述

产品名称	Anti-ILDR2抗体
描述	兔多克隆抗体to ILDR2
宿主	Rabbit
经测试应用	适用于: ICC/IF, IHC-P, WB
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to amino acids 362-493 of Human ILDR2 (Q71H61).
阳性对照	RT4 and U251 MG lysates; Human appendix tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121369** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment

应用	Ab评论	说明
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		1/250 - 1/500.

## 靶标


### 序列相似性

Belongs to the immunoglobulin superfamily. LISCH7 family.  
Contains 1 Ig-like V-type (immunoglobulin-like) domain.

### 细胞定位


Membrane.

## 图片




Immunofluorescent staining of Human cell line U-251MG shows positivity in nuclei but not nucleoli. Recommended concentration of ab121369 1-4  $\mu\text{g/ml}$ . Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-ILDR2 antibody (ab121369)



ab121369 at 1/75 dilution staining ILDR2 in Paraffin-embedded Human Appendix tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ILDR2 antibody (ab121369)



Western blot - Anti-ILDR2 antibody (ab121369)

**All lanes :** Anti-ILDR2 antibody (ab121369) at 1/250 dilution

**Lane 1 :** RT4 lysate

**Lane 2 :** U251 MG lysate

**Lane 3 :** Human plasma

**Lane 4 :** Human liver lysate

**Lane 5 :** Human tonsil lysate

Developed using the ECL technique.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors