

Product datasheet

KN-93 (water soluble) ab120980

1 References 1 图像

概述


| | |
|--------------|---|
| 产品名称 | KN-93 (water soluble) |
| 描述 | Potent, specific CaMK II抑制剂 |
| 生物学描述 | Potent, specific and cell-permeable calmodulin-dependent protein kinase (CaMK) II inhibitor (IC ₅₀ = 370 nM). Also a potent and direct extracellular open channel voltage-gated potassium channel blocker (IC ₅₀ = 307 nM for Kv1.5); independent of CaM kinase II inhibition. Water soluble. |
| 纯度 | > 99% |

性能

化学名称 *N*-[2-[[[(2*E*)-3-(4-Chlorophenyl)-2-propen-1-yl]methylamino]methyl]phenyl]-*N*-(2-hydroxyethyl)-4-methoxybenzenesulfonamide phosphate

分子量 599.03

化学结构



分子式 C₂₆H₂₉ClN₂O₄S.H₃PO₄

CAS编号 1188890-41-6

PubChem识别号 16760530

存放说明 Store at -20°C. It is important to note that this product is reported to be light sensitive. Store In the Dark. Store under desiccating conditions.

溶解度概述 Soluble in water to 100 mM and in DMSO to 100 mM

处理 Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour.

Need more advice on solubility, usage and handling? Please visit our [frequently asked questions \(FAQ\) page](#) for more details.

SMILES

Clc3ccc(/C=C/CN(C)Cc2ccccc2N(CCO)S(=O)(=O)c1ccc(OC)cc1)cc3.O=P(O)(O)O

来源

Synthetic


应用

Our [Abpromise guarantee](#) covers the use of **ab120980** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|--------------------|------|--|
| Functional Studies | | Use at an assay dependent concentration. |

图片



[ab40754](#) staining cyclin D1 in MCF7 cells treated with KN-93 (water soluble) (ab120980), by ICC/IF. Decrease in cyclin D1 expression correlates with increased concentration of KN-93 (water soluble), as described in literature.

The cells were incubated at 37°C for 24 hours in media containing different concentrations of ab120980 (KN-93 (water soluble)) in DMSO, fixed with 100% methanol for 5 minutes at -20°C and blocked with PBS containing 10% goat serum, 0.3 M glycine, 1% BSA and 0.1% tween for 2h at room temperature. Staining of the treated cells with [ab40754](#) (1/100 dilution) was performed overnight at 4°C in PBS containing 1% BSA and 0.1% tween. A DyLight 488 goat anti-rabbit polyclonal antibody ([ab96899](#)) at 1/250 dilution was used as the secondary antibody. Nuclei were counterstained with DAPI and are shown in blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE, NOT FOR USE IN HUMANS"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team